

# Product Datasheet

## ImmPRESS(TM) Excel Amplified anti-Rabbit IgG HRP Polymer Staining Kit (Colorimetric) MP-7601-NB-15ml

Unit Size: 15 ml

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/MP-7601-NB](http://www.novusbio.com/MP-7601-NB)

Updated 11/13/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/MP-7601-NB](http://www.novusbio.com/reviews/destination/MP-7601-NB)



**MP-7601-NB-15ml****ImmPRESS(TM) Excel Amplified anti-Rabbit IgG HRP Polymer Staining Kit (Colorimetric)**


<b>Product Information</b>	
<b>Unit Size</b>	15 ml
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Store at 4C. Do not freeze.
<b>Conjugate</b>	HRP Polymer

<b>Product Description</b>	
<b>Description</b>	<b>Features:</b> Complete staining kits Highest sensitivity Precise, sharp antigen localization Low background Amplified polymer technology Convenient prediluted, ready-to-use reagents
<b>Species</b>	Rabbit
<b>Specificity/Sensitivity</b>	ImmPRESS(TM) Excel Staining Kit (15 ml) will stain approximately 75-150 sections and an ImmPRESS(TM) Excel Staining Kit (50 ml) will stain approximately 250-500 sections.
<b>Kit Components</b>	BLOXALL(TM) Endogenous HRP/AP Blocking Solution (r.t.u.), 2.5% Normal serum for blocking non-specific binding (r.t.u.), Amplifier(TM) Antibody (either goat anti-mouse IgG or goat anti-rabbit IgG), ImmPRESS(TM) Excel Amplified HRP Polymer Reagent (anti-goat IgG), ImmPACT(TM) DAB EqV Substrate
<b>Notes</b>	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.

<b>Product Application Details</b>	
<b>Applications</b>	ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, In-situ Hybridization
<b>Recommended Dilutions</b>	ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, In-situ Hybridization

**Images**

Immunohistochemistry: ImmPRESS(TM) Excel anti-Rabbit IgG Staining Kit (Colorimetric) [MP-7601-NB]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/MP-7601-NB](http://www.novusbio.com/reviews/submit/MP-7601-NB)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

