

Product Datasheet

ImmPRESS(TM) Duet anti-Mouse IgG, anti-Rabbit IgG HRP/AP Polymer Staining Kit (Colorimetric) MP-7724-NB

Unit Size: 1 Kit

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MP-7724-NB

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MP-7724-NB



MP-7724-NB

ImmPRESS(TM) Duet anti-Mouse IgG, anti-Rabbit IgG HRP/AP Polymer Staining Kit (Colorimetric)

| Product Information | |
|-----------------------------|--|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Store at 4C. Do not freeze. |
| Conjugate | HRP/AP Polymer |
| Product Description | |
| Description | Features: Complete Kits for IHC Double Labeling Convenient Dual Specificity (anti-mouseIgG/rabbit IgG) ImmPRESS™ HRP/AP Polymer Reagent High sensitivity and Specificity Non-Biotin Based Reagents |
| Species | Mouse, Rabbit |
| Specificity/Sensitivity | The ImmPRESS(TM) Duet Double Staining Polymer Kit is intended for use on non-rodent tissue specimens. Reliable and reproducible results are obtained on tissue sections where the two target antigens do not over-lap (co-localize) in the same structure of the same cell, but rather are expressed in different cell compartments or different cell types. Also included, BLOXALL(TM) Blocking Solution abolishes all endogenous peroxidase and alkaline phosphatase activity in a single 10 minute incubation. In addition, 2.5% Normal Serum blocks all non-specific binding of detection reagents to ensure crisp, specific staining. All reagents are supplied in stabilized, ready-to-use (R.T.U.) solutions for maximum convenience, and consistency between assays. The ImmPRESS Duet Double Staining Polymer Kits offer a complete and simple solution for visualizing target antigens in cell and tissue preparations and contains sufficient reagents to stain approximately 150 sections. |
| Kit Components | 15 ml BLOXALL™ Blocking Solution, 15 ml 2.5% Normal Serum, 15 ml ImmPRESS™ Duet HRP/AP Polymer Reagent (mixture of ImmPRESS™ HRP Anti-Mouse IgG and ImmPRESS™ AP Anti-Rabbit IgG), ImmPACT™ DAB EqV Substrate (HRP, brown), ImmPACT™ Vector® Red Substrate (AP, magenta) |
| Notes | This product is manufactured by Vector Laboratories and distributed by Novus Biologicals. |
| Product Application Details | |
| Applications | Immunohistochemistry, Immunocytochemistry |
| Recommended Dilutions | Immunohistochemistry, Immunocytochemistry |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MP-7724-NB

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

