

Product Datasheet

Fibronectin Antibody NB110-41562-0.1ml

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-41562

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-41562



NB110-41562-0.1ml**Fibronectin Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.4), 0.2% BSA, Tween-20


Product Description	
Description	Novus Biologicals Rabbit Fibronectin Antibody (NB110-41562) is a polyclonal antibody validated for use in IHC. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2335
Gene Symbol	FN1
Species	Human
Reactivity Notes	Human.
Marker	Mesenchymal Cells Marker
Specificity/Sensitivity	This reacts with human fibronectin in liver, tonsil, skin and kidney. Traces of contaminating antibodies have been removed by solid-phase absorption. Only the fibronectin arch was observed in crossed immunoelectrophoresis.
Immunogen	Fibronectin isolated from a pool of normal human plasma.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Immunohistochemistry 1:50-1:100, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen 1:50-1:250
Application Notes	IHC-P: recommended pretreatment of EDTA buffer, pH 8.0. Recommended incubation time of 30 min at RT.




Images

Immunohistochemistry-Paraffin: Fibronectin Antibody [NB110-41562] - Formalin fixed paraffin embedded human kidney stained with Fibronectin antibody.(NB110-41562)



Immunohistochemistry-Paraffin: Fibronectin Antibody [NB110-41562] - Formalin fixed paraffin embedded human kidney stained with Fibronectin antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-41562-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-91258PEP	Fibronectin Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-41562

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

