

# Product Datasheet

## MYH6 Antibody (3-48) [DyLight 550] NB300-284R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB300-284R](http://www.novusbio.com/NB300-284R)

Updated 12/30/2019 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB300-284R](http://www.novusbio.com/reviews/destination/NB300-284R)



**NB300-284R**

MYH6 Antibody (3-48) [DyLight 550]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                             |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Clone                       | 3-48                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Isotype                     | IgG1 Kappa                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Conjugate                   | DyLight 550                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Purity                      | Protein A or G purified                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Gene ID                     | 4624                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Gene Symbol                 | MYH6                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Species                     | Human, Mouse, Rat, Canine                                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactivity Notes            | Canine reactivity reported in (PMID: 23577316). Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. Predicted to react with Rabbit, Bovine. |
| Specificity/Sensitivity     | This antibody recognizes alpha and beta type heavy chains of light meromyosin fragments (LMM) of cardiac myosin heavy chains. It also reacts with LMM fragments of cardiac myosin and beta heavy chain of slow human skeletal muscle. The antibody does not crossreact with human myosin light chains.                                                                                                                                  |
| Immunogen                   | Human ventricular myosin                                                                                                                                                                                                                                                                                                                                                                                                                |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                                       |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                               |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB300-284R](http://www.novusbio.com/reviews/submit/NB300-284R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

