

# Product Datasheet

## u-Plasminogen Activator/Urokinase Antibody (U-6) [FITC] NBP2-23504F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-23504F](http://www.novusbio.com/NBP2-23504F)

Updated 9/10/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-23504F](http://www.novusbio.com/reviews/destination/NBP2-23504F)



**NBP2-23504F**

u-Plasminogen Activator/Urokinase Antibody (U-6) [FITC]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | U-6                                                                                         |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | FITC                                                                                        |
| Purity              | Protein A or G purified                                                                     |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                                                                                                                                       |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                                                                                                                                                                                                                                                                 |
| Gene ID                 | 5328                                                                                                                                                                                                                                                                                                                                                                                                  |
| Gene Symbol             | PLAU                                                                                                                                                                                                                                                                                                                                                                                                  |
| Species                 | Human                                                                                                                                                                                                                                                                                                                                                                                                 |
| Reactivity Notes        | Human                                                                                                                                                                                                                                                                                                                                                                                                 |
| Specificity/Sensitivity | Specific for human uPA. No reaction with human tissue plasminogen activator (tPA) is seen when tested by ELISA, immunoblotting and enzyme inhibition. No reaction with any other human plasma proteins is seen when tested by immunoblotting. Binds single and two-chain uPA, uPA/PAI-1 complex and uPA/uPAR complex. A reaction is seen with the amino terminal fragment (ATF) but not with LMW uPA. |
| Immunogen               | Native human uPA                                                                                                                                                                                                                                                                                                                                                                                      |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | ELISA, Immunohistochemistry, Immunohistochemistry-Frozen               |
| Recommended Dilutions       | ELISA, Immunohistochemistry, Immunohistochemistry-Frozen               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-23504F](http://www.novusbio.com/reviews/submit/NBP2-23504F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

