

Product Datasheet

Pancreatic trypsin inhibitor Antibody (35) - BSA Free NBP2-23590-0.1ml

Unit Size: 0.1ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-23590

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-23590



NBP2-23590-0.1ml

Pancreatic trypsin inhibitor Antibody (35) - BSA Free

Product Information	
Unit Size	0.1ml
Concentration	1 mg/ml
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	35
Preservative	0.09% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	10mM Phosphate (pH 7.4) and 0.5M NaCl

Product Description	
Description	Novus Biologicals Mouse Pancreatic trypsin inhibitor Antibody (35) - BSA Free (NBP2-23590) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	404172
Gene Symbol	PTI
Species	Bovine
Specificity/Sensitivity	Clone 35 binds bovine aprotinin immobilized on a solid phase or free in solution.
Immunogen	Aprotinin from bovine lung

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:100-1:2000, ELISA 1:100-1:2000
Application Notes	ELISA: NBP2-23590 (as capture antibody coated at 1 ug/mL) forms an antibody pair with NBP2-23588 (as biotinylated detection antibody) for sandwich ELISA of aprotinin at sample concentrations of 4-250 ug/mL. Both urine and plasma samples diluted 1/1000 can be assayed against standards in buffer. PBS-based buffers are unsuitable for this assay and 0.2 M glycine sodium buffer, pH 8.8, is recommended for washing and dilution buffers.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-23590-0.1ml

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP3-30997	Bovine Pancreatic trypsin inhibitor ELISA Kit (Colorimetric)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-23590

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

