

# Product Datasheet

## epithelial Sodium Channel gamma Antibody (3c7) [DyLight 488]

**NBP2-41371G**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-41371G](http://www.novusbio.com/NBP2-41371G)

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-41371G](http://www.novusbio.com/reviews/destination/NBP2-41371G)



**NBP2-41371G**

epithelial Sodium Channel gamma Antibody (3c7) [DyLight 488]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 3c7                                                                                         |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1 Kappa                                                                                  |
| Conjugate           | DyLight 488                                                                                 |
| Purity              | Protein A or G purified                                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                                                                               |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                         |
| Gene ID                 | 6340                                                                                                                                                          |
| Gene Symbol             | SCNN1G                                                                                                                                                        |
| Species                 | Human                                                                                                                                                         |
| Specificity/Sensitivity | Antibody is specific for the inhibitory tract of human gamma-ENaC subunit.                                                                                    |
| Immunogen               | The inhibitory peptide from the human gamma-ENaC subunit.<br>EAESWNSVSEGKQPRFSHRIPLC corresponding to amino acid residue 139-160 of human gamma-ENaC subunit. |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                             |

| Product Application Details |                                                                                                                                               |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Sandwich ELISA, Sandwich ELISA Capture |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Sandwich ELISA, Sandwich ELISA Capture |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                        |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-41371G**

---

|                    |                                                                               |
|--------------------|-------------------------------------------------------------------------------|
| NBP1-43319G-0.5ml  | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 488]                   |
| NBP2-56532PEP      | epithelial Sodium Channel gamma Recombinant Protein Antigen                   |
| H00023327-P01-10ug | Recombinant Human NEDD4L GST (N-Term) Protein                                 |
| NBP3-31076         | Human epithelial Sodium Channel gamma - Ready-To-Use ELISA Kit (Colorimetric) |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-41371G](http://www.novusbio.com/reviews/submit/NBP2-41371G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

