

# Product Datasheet

## **Triosephosphate isomerase Antibody (1111) [DyLight 488] NBP2-43824G**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-43824G](http://www.novusbio.com/NBP2-43824G)

Updated 8/1/2019 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-43824G](http://www.novusbio.com/reviews/destination/NBP2-43824G)



**NBP2-43824G**

Triosephosphate isomerase Antibody (1111) [DyLight 488]

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                                         |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | 1111                                                                                                             |
| Preservative        | 0.05% Sodium Azide                                                                                               |
| Isotype             | IgG1                                                                                                             |
| Conjugate           | DyLight 488                                                                                                      |
| Purity              | Protein G purified                                                                                               |
| Buffer              | 50mM Sodium Borate                                                                                               |

| Product Description |                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                                                                                    |
| Gene ID             | 7167                                                                                                                                                                     |
| Gene Symbol         | TPI1                                                                                                                                                                     |
| Species             | Human, Mouse, Rat                                                                                                                                                        |
| Immunogen           | Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human Triosephosphate isomerase. The exact sequence is proprietary. |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
[canada.inquires@bio-techne.com](mailto:canada.inquires@bio-techne.com)

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
[info.EMEA@bio-techne.com](mailto:info.EMEA@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [technical@novusbio.com](mailto:technical@novusbio.com)  
Orders: [orders@novusbio.com](mailto:orders@novusbio.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-43824G](http://www.novusbio.com/reviews/submit/NBP2-43824G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

