

# Product Datasheet

## EGFR Antibody (Matuzumab) [PE/Cy5.5] - Chimeric NBP2-52685PECY55

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-52685PECY55](http://www.novusbio.com/NBP2-52685PECY55)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-52685PECY55](http://www.novusbio.com/reviews/destination/NBP2-52685PECY55)



**NBP2-52685PECY55**

EGFR Antibody (Matuzumab) [PE/Cy5.5] - Chimeric

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark. Do not freeze.                                                     |
| Clonality           | Monoclonal                                                                                  |
| Clone               | Matuzumab                                                                                   |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG Kappa                                                                                   |
| Conjugate           | PE/Cy5.5                                                                                    |
| Purity              | Protein A purified                                                                          |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene ID                 | 1956                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Gene Symbol             | EGFR                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Species                 | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactivity Notes        | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Specificity/Sensitivity | This antibody binds the extracellular domain III of the human EGFR and does not cross react with murine EGFR. It was also reported that residues S460/G461 in EGFR domain III are essential components of the epitope. This antibody binds at a site distinct from the EGF binding site, preventing conformational re-arrangement required for dimerization. This antibody can only recognize the active form of EGFR. The antibody binds to both glyco- and aglyco receptor forms. |
| Immunogen               | The parental mouse antibody was generated by immunizing BALB/c mice intraperitoneally with A431 cells in phosphate buffered saline (PBS). Later on the humanized version of the antibody was created by grafting CDRs of the murine antibody onto human constant regions.                                                                                                                                                                                                           |

| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Flow Cytometry                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Recommended Dilutions       | Flow Cytometry                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at <a href="mailto:technical@novusbio.com">technical@novusbio.com</a> if you have any questions. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-52685PECY55**

---

|               |                                  |
|---------------|----------------------------------|
| NBP2-52952PEP | EGFR Recombinant Protein Antigen |
| 236-EG-200    | EGF [Unconjugated]               |
| DYC1095B-2    | EGFR [Biotin]                    |
| DVE00         | VEGF [HRP]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-52685PECY55](http://www.novusbio.com/reviews/submit/NBP2-52685PECY55)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

