

Product Datasheet

beta-III Tubulin Antibody (CL5813) - BSA Free NBP2-61430

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61430

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61430



NBP2-61430


beta-III Tubulin Antibody (CL5813) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL5813
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Product Description	
Description	Novus Biologicals Mouse beta-III Tubulin Antibody (CL5813) - BSA Free (NBP2-61430) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10381
Gene Symbol	TUBB3
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
Immunogen	This antibody was developed using a synthetic peptide derived from Q13509, with the exact immunogen sequence remaining proprietary.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:50000 - 1:100000, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:50000 - 1:100000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation Permeabilization: PFA/Triton X-100




Images

Western Blot: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Analysis in human cerebral cortex tissue.




Immunocytochemistry/Immunofluorescence: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of A549 cells showing specific staining on microtubules in green. Nuclear probes are visualized in blue. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of rat cerebral cortex shows strong positivity in neuronal cell bodies and in neuropil.




Western Blot: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Analysis in human cell line SH-SY5Y.




Immunocytochemistry/Immunofluorescence: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of SH-SY5Y cells using the Anti-TUBB3 monoclonal antibody, showing specific staining on microtubules in green. Nuclear probe is visualized in blue.




Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of human cerebral cortex shows strong positivity in neuropil as well as in neurons.




Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of human cerebellum shows strong positivity in both neurons and neuropil.




Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of human prostate shows strong positivity in the peripheral nerves.



Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of human skeletal muscle shows absence of positivity in striated muscle fibers as expected (negative control).



Immunohistochemistry-Paraffin: beta-III Tubulin Antibody (CL5813) [NBP2-61430] - Staining of mouse embryo E11 shows strong positivity in the developing central and peripheral nervous system.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61430

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-22901	Recombinant Human beta-III Tubulin His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61430

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

