

Product Datasheet

P2Y1/P2RY1 Antibody (4G5B7) - BSA Free NBP2-61671

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61671

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61671



NBP2-61671

P2Y1/P2RY1 Antibody (4G5B7) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4G5B7
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	42 kDa

Product Description	
Description	Novus Biologicals Mouse P2Y1/P2RY1 Antibody (4G5B7) - BSA Free (NBP2-61671) is a monoclonal antibody validated for use in IHC, WB, ELISA and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5028
Gene Symbol	P2RY1
Species	Human
Immunogen	Purified recombinant fragment of human P2Y1/P2RY1 (AA: extra mix) expressed in E. Coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Western Blot 1:100-1:2000, Flow Cytometry 1:10-1:1000, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:1000, Immunohistochemistry-Paraffin 1:200 - 1:1000

Images


Western Blot: P2Y1/P2RY1 Antibody (4G5B7) [NBP2-61671] - Analysis using P2RY1 mAb against human P2RY1 (AA: extra mix) recombinant protein. (Expected MW is 13 kDa)




Immunohistochemistry-Paraffin: P2Y1/P2RY1 Antibody (4G5B7) [NBP2-61671] - Analysis of rectum cancer tissues using P2RY1 mouse mAb with DAB staining.




Flow Cytometry: P2Y1/P2RY1 Antibody (4G5B7) [NBP2-61671] - Analysis of Hela cells using P2RY1 mouse mAb (green) and negative control (red).




Western Blot: P2Y1/P2RY1 Antibody (4G5B7) [NBP2-61671] - Analysis using P2RY1 mAb against HEK293 (1) and P2RY1 (AA: extra mix)-hlgGfc transfected HEK293 (2) cell lysate.



Immunohistochemistry-Paraffin: P2Y1/P2RY1 Antibody (4G5B7) [NBP2-61671] - Analysis of bladder cancer tissues using P2RY1 mouse mAb with DAB staining.



ELISA: P2Y1/P2RY1 Antibody (4G5B7) [NBP2-61671] - Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61671

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
H00005028-Q01-10ug	Recombinant Human P2Y1/P2RY1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61671

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

