

Product Datasheet

Histone H3 [p Thr45] Antibody - BSA Free NBP2-62795

Unit Size: 50 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-62795

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-62795



NBP2-62795

Histone H3 [p Thr45] Antibody - BSA Free

Product Information

Unit Size	50 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide and 0.05% ProClin 300
Isotype	IgG
Purity	Affinity purified
Buffer	PBS
Target Molecular Weight	15 kDa

Product Description

Description	Novus Biologicals Rabbit Histone H3 [p Thr45] Antibody - BSA Free (NBP2-62795) is a polyclonal antibody validated for use in ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	126961
Gene Symbol	H3C14
Species	Human
Immunogen	This Histone H3 [p Thr45] antibody was raised against histone H3 containing the phosphorylated threonine 45 (H3T45p)


Product Application Details

Applications	Dot Blot, ELISA
Recommended Dilutions	ELISA 1:5000 - 1:10000, Dot Blot 1:20000




Images

ELISA: Histone H3 [p Thr45] Antibody [NBP2-62795] - To determine the titer of the antibody, analysis was performed using a serial dilution of the Diagenode antibody directed against H3T45p in antigen coated wells. The antigen used was a peptide containing the histone modification of interest. By plotting the absorbance against the antibody dilution, the titer of the antibody was estimated to be 1:350,000.



Dot Blot: Histone H3 [p Thr45] Antibody [NBP2-62795] - To test the cross reactivity of the antibody against H3T45p, an analysis was performed with peptides containing other histone modifications and the unmodified H3T45. One hundred to 0.2 pmol of the respective peptides were spotted on a membrane. This figure shows a high specificity of the antibody for the modification of interest.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-62795

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB21-1162PEP	Histone H3 [Dimethyl Lys23] Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-62795

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

