

Product Datasheet

Histone H3 [p Ser10] Antibody (SA31-01) NBP2-67305

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-67305

Updated 2/9/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-67305



NBP2-67305

Histone H3 [p Ser10] Antibody (SA31-01)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	SA31-01
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	TBS (pH7.4), 0.05% BSA, 40% Glycerol
Target Molecular Weight	15 kDa

Product Description	
Host	Rabbit
Gene Symbol	H3C14
Species	Human, Mouse, Rat
Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser10 of human Histone H3. (SwissProt: P68431 Human; SwissProt: P84243 Human; SwissProt: Q16695 Human; SwissProt: Q6NXT2 Human; SwissProt: Q71DI3 Human; SwissProt: P84244 Mouse; SwissProt: P84245 Rat)

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:5000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunoprecipitation, Immunohistochemistry-Paraffin 1:50-1:200

Images


Western Blot: Histone H3 [p Ser10] Antibody (SA31-01) [NBP2-67305] - Analysis of Histone H3 (p Ser10) on different lysates with Rabbit antiHistone H3 (p Ser10) antibody at 1/500 dilution. Lane 1: K562 cell lysate Lane 2: U937 cell lysate Lysates/proteins at 10 ug/Lane. Predicted band size: 15 kDa Observed band size: 15 kDa Exposure time: 2 minutes; 15% SDS-PAGE gel. Proteins were transferred to a PVDF membrane and blocked with 5% NFDm/TBST for 1 hour at room temperature. The primary antibody at 1/500 dilution was used in 5% NFDm/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:300,000 dilution was used for 1 hour at room temperature.




Immunocytochemistry/Immunofluorescence: Histone H3 [p Ser10] Antibody (SA31-01) [NBP2-67305] - Staining Phospho-Histone H3(S10) in Ags cells (green). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: Histone H3 [p Ser10] Antibody (SA31-01) [NBP2-67305] - Analysis of paraffin-embedded mouse testis tissue using anti-Phospho-Histone H3(S10) antibody. Counter stained with hematoxylin.




Western Blot: Histone H3 [p Ser10] Antibody (SA31-01) [NBP2-67305] - Analysis of Phospho-Histone H3(S10) on HeLa cell lysates using anti-Phospho-Histone H3(S10) antibody at 1/1,000 dilution. Positive control: Lane 1: Untreated HeLa cell lysate Lane 2: HeLa cell lysate treated with calyculin A Lane 3: HeLa cell lysate treat



Immunocytochemistry/Immunofluorescence: Histone H3 [p Ser10] Antibody (SA31-01) [NBP2-67305] - Staining Phospho-Histone H3(S10) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunocytochemistry/Immunofluorescence: Histone H3 [p Ser10]
Antibody (SA31-01) [NBP2-67305] - Staining Phospho-Histone H3(S10)
in HepG2 cells (green). Cells were fixed in paraformaldehyde,
permeabilised with 0.25% Triton X100/PBS.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-67305

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

