

Product Datasheet

PIM1 Antibody (ST0513) NBP2-67528

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 7

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-67528

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-67528



NBP2-67528

PIM1 Antibody (ST0513)

Product Information

Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	ST0513
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	TBS (pH7.4), 0.05% BSA, 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit PIM1 Antibody (ST0513) (NBP2-67528) is a recombinant monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. Anti-PIM1 Antibody: Cited in 7 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5292
Gene Symbol	PIM1
Species	Human, Mouse
Immunogen	Synthetic peptide within human PIM1 aa 20-60. (SwissProt: P11309 Human; SwissProt: P06803 Mouse)

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Western Blot 1:500-1:2000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin 1:50-1:400, Immunohistochemistry-Frozen
Application Notes	Use in IHC-Fr reported in scientific literature (PMID:35896115).

Images


Western Blot: PIM1 Antibody (ST0513) [NBP2-67528] - Analysis of PIM1 on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:5,000 dilution was used for 1 hour at room temperature. Positive control: Lane 1: BT-20 cell lysate Lane 2: HepG2 cell lysate Lane 3: SW480 cell lysate




Immunocytochemistry/Immunofluorescence: PIM1 Antibody (ST0513) [NBP2-67528] - Staining PIM1 in SW480 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: PIM1 Antibody (ST0513) [NBP2-67528] - Analysis of human tonsil tissue using anti-PIM1 antibody. Counter stained with hematoxylin.




Immunocytochemistry/Immunofluorescence: PIM1 Antibody (ST0513) [NBP2-67528] - Staining PIM1 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunocytochemistry/Immunofluorescence: PIM1 Antibody (ST0513) [NBP2-67528] - Staining PIM1 in BT-20 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: PIM1 Antibody (ST0513) [NBP2-67528]
- Analysis of human breast carcinoma tissue using anti-PIM1 antibody.
Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: PIM1 Antibody (ST0513) [NBP2-67528]
- Analysis of human colon cancer tissue using anti-PIM1 antibody.
Counter stained with hematoxylin.



Immunohistochemistry-Paraffin: PIM1 Antibody (ST0513) [NBP2-67528]
- Analysis of human gastric carcinoma tissue using anti-PIM1 antibody.
Counter stained with hematoxylin.



Immunohistochemistry-Paraffin: PIM1 Antibody (ST0513) [NBP2-67528]
- Analysis of human spleen tissue using anti-PIM1 antibody. Counter
stained with hematoxylin.



Publications

Zhu L, Li H, Peng X et Al. Beneficial mechanisms of dimethyl fumarate in autoimmune uveitis: insights from single-cell RNA sequencing J Neuroinflammation 2024-04-29 [PMID: 38684986]

Miruna Cristian, Mariana Așchie, Anca-Florentina Mitroi, Mariana Deacu, Mădălina Boșoteanu, Gabriela-Izabela Bălătescu, Andreea-Georgiana Stoica, Anca-Antonela Nicolau, Manuela Enciu, Ana-Maria Crețu, Andreea-Daniela Caloian, Cristian-Ionuț Orășanu, Ionuț Poinăreanu The impact of MYD88 and PIM1 in mature large B-cell non-Hodgkin lymphomas: Defining element of their evolution and prognosis Medicine 2024-02-09 [PMID: 38335426]

Li Y, Shi L, Zhao F et al. PIM1 attenuates cisplatin-induced AKI by inhibiting Drp1 activation Cellular signalling 2024-01-01 [PMID: 37967691] (ICC/IF, IHC, Mouse)

Luo Y, Zhu J, Zhao F et al. PIM1 attenuates renal ischemia-reperfusion injury by inhibiting ASK1-JNK/P38 International immunopharmacology 2022-12-10 [PMID: 36513021] (IHC-P, Mouse)

Li H, Xie L, Zhu L et al. Multicellular immune dynamics implicate PIM1 as a potential therapeutic target for uveitis Nature communications 2022-10-04 [PMID: 36195600] (FLOW, IF/IHC, Mouse)

Li Z, Wang M, Tan J et al. Single-cell RNA sequencing depicts the local cell landscape in thyroid-associated ophthalmopathy Cell reports. Medicine 2022-07-19 [PMID: 35896115] (IHC-Fr, Human)

Szydlowski M, Debek S, Prochorec-Sobieszek M et al. PIM Kinases Promote Survival and Immune Escape in Primary Mediastinal Large B-Cell Lymphoma through Modulation of JAK-STAT and NF-kappaB Activity The American journal of pathology 2020-12-08 [PMID: 33307035]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-67528

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-67528

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

