

# Product Datasheet

## CAMLG Antibody (OTI1A3) [PerCP] NBP2-70331PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-70331PCP](http://www.novusbio.com/NBP2-70331PCP)

Updated 11/9/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-70331PCP](http://www.novusbio.com/reviews/destination/NBP2-70331PCP)



**NBP2-70331PCP**

CAMLG Antibody (OT11A3) [PerCP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | OT11A3                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | PerCP                                                                                       |
| Purity              | Immunogen affinity purified                                                                 |
| Buffer              | PBS                                                                                         |

| Product Description |                                                                                       |
|---------------------|---------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                 |
| Gene ID             | 819                                                                                   |
| Gene Symbol         | CAMLG                                                                                 |
| Species             | Human, Mouse, Rat                                                                     |
| Immunogen           | Full length recombinant protein of human CAMLG (NP_001736) produced in HEK293T cells. |

| Product Application Details |                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                              |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-70331PCP**

---

|               |                                          |
|---------------|------------------------------------------|
| NBP1-97005PCP | Mouse IgG1 Isotype Control (MG1) [PerCP] |
| NBP2-22725    | Recombinant Human CAMLG His Protein      |
| 210-TA-005    | TNF-alpha [Unconjugated]                 |
| NBP2-08179    | CAMLG Overexpression Lysate              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-70331PCP](http://www.novusbio.com/reviews/submit/NBP2-70331PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

