

# Product Datasheet

## **FGFR2 Antibody (OTI5C5) [Janelia Fluor® 646] NBP2-70731JF646**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-70731JF646](http://www.novusbio.com/NBP2-70731JF646)

Updated 11/9/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-70731JF646](http://www.novusbio.com/reviews/destination/NBP2-70731JF646)



**NBP2-70731JF646**

FGFR2 Antibody (OTI5C5) [Janelia Fluor® 646]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | OTI5C5                                                                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                                                                                                           |
| Conjugate                   | Janelia Fluor 646                                                                                                                                                                                                                                                                                                                               |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                                                     |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 2263                                                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | FGFR2                                                                                                                                                                                                                                                                                                                                           |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                               |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen                   | Full length recombinant protein of human FGF R2 (NP_000132) produced in HEK293T cells.                                                                                                                                                                                                                                                          |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                                                                                                                                                                                           |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                             |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-70731JF646**

---

|                 |                                          |
|-----------------|------------------------------------------|
| NBP2-59905-50ug | Recombinant Human FGFR2 hlgG-His Protein |
| 236-EG-200      | EGF [Unconjugated]                       |
| NBP2-04392      | FGFR2 Overexpression Lysate              |
| DVE00           | VEGF [HRP]                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-70731JF646](http://www.novusbio.com/reviews/submit/NBP2-70731JF646)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

