

Product Datasheet

IL-6 Antibody (OTI3G9) - Azide and BSA Free NBP2-71027

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-71027

Updated 5/2/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-71027



NBP2-71027

IL-6 Antibody (OTI3G9) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI3G9
Preservative	No Preservative
Reconstitution Instructions	Reconstitute with PBS (pH 7.3). To use this carrier-free antibody for conjugation experiments, another round of desalting is highly recommended.
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose
Target Molecular Weight	23.5 kDa

Product Description	
Host	Mouse
Gene ID	3569
Gene Symbol	IL6
Species	Human
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-212 of human IL-6 (NP_000591) produced in E.coli.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:50000, Immunohistochemistry 1:50, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin

Images


Western Blot: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - IL6 antibody (3G9) at 1:500000 + recombinant human IL6 antibody (2G3) at 1:2000




Immunocytochemistry/Immunofluorescence: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of paraffin-embedded Human endometrium tissue using anti-IL6 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of paraffin-embedded Carcinoma of Human pancreas tissue using anti-IL6 mouse monoclonal antibody.




Western Blot: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IL-6 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IL-6. Positive lysates LY424617 (100ug) and LC424617 (20ug).




Immunocytochemistry/Immunofluorescence: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Anti-IL-6 mouse monoclonal antibody immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IL-6.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human breast tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Adenocarcinoma of Human breast tissue using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human colon tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Adenocarcinoma of Human colon tissue using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human Kidney tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Carcinoma of Human kidney tissue using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human lung tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Carcinoma of Human lung tissue using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human Ovary tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.




Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Adenocarcinoma of Human ovary tissue using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.



Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human pancreas tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.



Immunohistochemistry-Paraffin: IL-6 Antibody (OTI3G9) - Azide and BSA Free [NBP2-71027] - Staining of Human prostate tissue within the normal limits using anti-IL-6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100 degrees C for 10min.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-71027

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
NBP2-35134-10ug	Recombinant Mouse IL-6 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-71027

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

