

# Product Datasheet

## KCTD14 Antibody (OTI1D11) [DyLight 550] NBP2-72173R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-72173R](http://www.novusbio.com/NBP2-72173R)

Updated 11/9/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-72173R](http://www.novusbio.com/reviews/destination/NBP2-72173R)



**NBP2-72173R**

KCTD14 Antibody (OTI1D11) [DyLight 550]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | OTI1D11                                                                                     |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | DyLight 550                                                                                 |
| Purity              | Immunogen affinity purified                                                                 |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                        |
|---------------------|----------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                  |
| Gene ID             | 65987                                                                                  |
| Gene Symbol         | KCTD14                                                                                 |
| Species             | Human                                                                                  |
| Immunogen           | Full length recombinant protein of human KCTD14 (NP_076419) produced in HEK293T cells. |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.      |

| Product Application Details |                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                              |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-72173R**

|               |                                                |
|---------------|------------------------------------------------|
| NBP1-97005R   | Mouse IgG1 Isotype Control (MG1) [DyLight 550] |
| NBP1-92043PEP | KCTD14 Recombinant Protein Antigen             |
| NBL1-12211    | KCTD14 Overexpression Lysate                   |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-72173R](http://www.novusbio.com/reviews/submit/NBP2-72173R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

