

Product Datasheet

Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free NBP2-72956

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-72956

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-72956




NBP2-72956**Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free**

Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI3B5
Preservative	No Preservative
Reconstitution Instructions	we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose
Target Molecular Weight	23.077 kDa
Product Description	
Description	Novus Biologicals Mouse Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free (NBP1-47878) is a monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	50674
Gene Symbol	NEUROG3
Species	Human
Immunogen	Full-length protein expressed in 293T cell transfected with human NEUROG3 expression vector
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry, Immunohistochemistry 1:50 - 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunoprecipitation, Immunohistochemistry-Paraffin




Images


Western Blot: Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free [NBP2-72956] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY Neurogenin-3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-Neurogeni




Immunohistochemistry: Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free [NBP2-72956] - Staining of paraffin-embedded Human lymph node tissue using anti-Neurogenin-3 mouse monoclonal antibody.




Flow Cytometry: Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free [NBP2-72956] - Neurogenin3 Staining of HeLa cells transiently transfected by pCMV6-ENTRY Neurogenin-3.



Immunohistochemistry: Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free [NBP2-72956] - Staining of paraffin-embedded Human liver tissue within the normal limits using anti-NEUROG3 mouse monoclonal antibody at 1:150 dilution.



Immunoprecipitation: Neurogenin-3 Antibody (OTI3B5) - Azide and BSA Free [NBP2-72956] - (Negative control: IP without adding anti-NEUROG3 antibody). For each experiment, 500ul of DDK tagged NEUROG3 overexpression lysates (at 1:5 dilution with HEK293T lysate), 2ug of anti-NEUROG3 antibody and 20ul (0.1mg) of goat anti-mouse conjugated magnetic beads were mixed and incubated overnight. After extensive wash to remove any non-specific binding, the immunoprecipitated products were analyzed with rabbit anti-DDK polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-72956

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
NBP2-51514-0.02mg	Recombinant Human Neurogenin-3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-72956

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

