

# Product Datasheet

## PIH1D2 Antibody (OTI4A10) [Biotin] NBP2-73403B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-73403B](http://www.novusbio.com/NBP2-73403B)

Updated 11/9/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-73403B](http://www.novusbio.com/reviews/destination/NBP2-73403B)



**NBP2-73403B**

PIH1D2 Antibody (OTI4A10) [Biotin]

| Product Information         |                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                     |
| Storage                     | Store at 4C in the dark.                                                                                        |
| Clonality                   | Monoclonal                                                                                                      |
| Clone                       | OTI4A10                                                                                                         |
| Preservative                | 0.05% Sodium Azide                                                                                              |
| Isotype                     | IgG1                                                                                                            |
| Conjugate                   | Biotin                                                                                                          |
| Purity                      | Immunogen affinity purified                                                                                     |
| Buffer                      | PBS                                                                                                             |
| Product Description         |                                                                                                                 |
| Host                        | Mouse                                                                                                           |
| Gene ID                     | 120379                                                                                                          |
| Gene Symbol                 | PIH1D2                                                                                                          |
| Species                     | Human                                                                                                           |
| Immunogen                   | Recombinant protein fragment corresponding to amino acids 1-315 of human PIH1D2 (NP_620144) produced in E.coli. |
| Product Application Details |                                                                                                                 |
| Applications                | Western Blot, Flow Cytometry, Immunohistochemistry                                                              |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry                                                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-73403B**

---

|                   |                                           |
|-------------------|-------------------------------------------|
| NBP1-96888        | Mouse IgG1 Isotype Control (MG1) [Biotin] |
| NBP2-51879-0.02mg | Recombinant Human PIH1D2 His Protein      |
| NBL1-14413        | PIH1D2 Overexpression Lysate              |
| NBP2-94190-0.1ml  | TIMM8B Antibody - Azide and BSA Free      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-73403B](http://www.novusbio.com/reviews/submit/NBP2-73403B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

