

# Product Datasheet

## STK39 Antibody (OTI4E3) [PE] NBP2-74406PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-74406PE](http://www.novusbio.com/NBP2-74406PE)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-74406PE](http://www.novusbio.com/reviews/destination/NBP2-74406PE)



**NBP2-74406PE**

STK39 Antibody (OTI4E3) [PE]

| Product Information         |  |
|-----------------------------|--|
| Unit Size                   | 0.1 ml   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Monoclonal   |
| Clone                       | OTI4E3   |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG2b  |
| Conjugate                   | PE   |
| Purity                      | Immunogen affinity purified  |
| Buffer                      | PBS  |
| Product Description         |  |
| Host                        | Mouse  |
| Gene ID                     | 27347  |
| Gene Symbol                 | STK39  |
| Species                     | Human, Mouse, Rat  |
| Immunogen                   | Full length recombinant protein of human STK39 (NP_037365) produced in HEK293T cells.  |
| Product Application Details |  |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin, CyTOF-ready |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-74406PE**

---

|                    |  |
|--------------------|--|
| NBP2-27231PE       | Mouse IgG2b Isotype Control (MPC-11) [PE]    |
| H00027347-Q01-10ug | Recombinant Human STK39 GST (N-Term) Protein |
| NBL1-16554         | STK39 Overexpression Lysate                  |
| NB600-284          | WNK4 Antibody - BSA Free                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-74406PE](http://www.novusbio.com/reviews/submit/NBP2-74406PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

