

Product Datasheet

TRAIP Antibody (OTI2D4) - Azide and BSA Free NBP2-74580

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-74580

Updated 10/9/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-74580



NBP2-74580

TRAIP Antibody (OTI2D4) - Azide and BSA Free

Product Information

Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI2D4
Preservative	No Preservative
Reconstitution Instructions	Reconstitute with PBS (pH 7.3). To use this carrier-free antibody for conjugation experiments, another round of desalting is highly recommended.
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose
Target Molecular Weight	53.1 kDa

Product Description


Host	Mouse
Gene ID	10293
Gene Symbol	TRAIP
Species	Human
Immunogen	Human recombinant protein fragment corresponding to amino acids 67-313 of human TRAIP (NP_005870) produced in E.coli.

Product Application Details


Applications	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, CyTOF-ready

Images


Western Blot: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRAIP (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRAIP.




Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Human Kidney tissue using anti-TRAIP mouse monoclonal antibody.




Flow Cytometry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostaining by anti-TRAIP antibody, and then analyzed by flow cytometry.




Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-TRAIP mouse monoclonal antibody.




Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-TRAIP mouse monoclonal antibody.




Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-TRAIP mouse monoclonal antibody.




Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Carcinoma of Human kidney tissue using anti-TRAIP mouse monoclonal antibody.



Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Carcinoma of Human liver tissue using anti-TRAIP mouse monoclonal antibody.




Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Carcinoma of Human lung tissue using anti-TRAIP mouse monoclonal antibody.



Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Human colon tissue using anti-TRAIP mouse monoclonal antibody.



Immunohistochemistry: TRAIP Antibody (OTI2D4) - Azide and BSA Free [NBP2-74580] - Staining of paraffin-embedded Human liver tissue using anti-TRAIP mouse monoclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-74580

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NB100-56180PEP	TRAIP Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-74580

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

