

# Product Datasheet

## CLIP170 Antibody NBP2-75095

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-75095](http://www.novusbio.com/NBP2-75095)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-75095](http://www.novusbio.com/reviews/destination/NBP2-75095)



**NBP2-75095****CLIP170 Antibody****Product Information**

|                                |                                          |
|--------------------------------|------------------------------------------|
| <b>Unit Size</b>               | 100 ug                                   |
| <b>Concentration</b>           | 0.5 mg/ml                                |
| <b>Storage</b>                 | Store at -20C. Avoid freeze-thaw cycles. |
| <b>Clonality</b>               | Polyclonal                               |
| <b>Preservative</b>            | 0.02% Sodium Azide                       |
| <b>Isotype</b>                 | IgG                                      |
| <b>Purity</b>                  | Antigen Affinity-purified                |
| <b>Buffer</b>                  | TBS (pH 7.3), 0.5% BSA                   |
| <b>Target Molecular Weight</b> | 161 kDa                                  |

**Product Description**


|                    |                                                                                                                                                                                      |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Novus Biologicals Goat CLIP170 Antibody (NBP2-75095) is a polyclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>        | Goat                                                                                                                                                                                 |
| <b>Gene ID</b>     | 6249                                                                                                                                                                                 |
| <b>Gene Symbol</b> | CLIP1                                                                                                                                                                                |
| <b>Species</b>     | Human                                                                                                                                                                                |
| <b>Immunogen</b>   | Peptide with sequence C-PHSTHHGSRGEERP., from the C Terminus of the protein sequence according to NP_002947.1; NP_937883.1.                                                          |

**Product Application Details**

|                              |                                          |
|------------------------------|------------------------------------------|
| <b>Applications</b>          | Western Blot, Peptide ELISA              |
| <b>Recommended Dilutions</b> | Western Blot 1:100-1:2000, Peptide ELISA |

**Images**

Western Blot: CLIP170 Antibody [NBP2-75095] - (1ug/ml) biotinylated staining of Human Breast Cancer lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence, using streptavidin-HRP and using NAP blocker as a substitute for skimmed milk.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-75095**

|                 |                                                                        |
|-----------------|------------------------------------------------------------------------|
| HAF017          | Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| HAF109          | Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| NB410-28088-1mg | Goat IgG Isotype Control                                               |
| NBP2-61143PEP   | CLIP170 Recombinant Protein Antigen                                    |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-75095](http://www.novusbio.com/reviews/submit/NBP2-75095)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)



