

# Product Datasheet

## Human COL9A3 ELISA Kit (Colorimetric) NBP2-75887

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-75887](http://www.novusbio.com/NBP2-75887)

Updated 9/23/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-75887](http://www.novusbio.com/reviews/destination/NBP2-75887)



**NBP2-75887****Human COL9A3 ELISA Kit (Colorimetric)****Product Information**

|                      |                                                                                                          |
|----------------------|----------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 1 Kit                                                                                                    |
| <b>Concentration</b> | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| <b>Storage</b>       | Storage of components varies. See protocol for specific instructions.                                    |

**Product Description**


|                                |                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>                 | 1299                                                                                                                                                                                                                                                                                                                                                               |
| <b>Gene Symbol</b>             | COL9A3                                                                                                                                                                                                                                                                                                                                                             |
| <b>Species</b>                 | Human                                                                                                                                                                                                                                                                                                                                                              |
| <b>Specificity/Sensitivity</b> | This kit recognizes Human COL9a3 in samples. No significant cross-reactivity or interference between Human COL9a3 and analogues was observed.                                                                                                                                                                                                                      |
| <b>Kit Components</b>          | Micro ELISA Plate(Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100×), Concentrated HRP Conjugate (100×), Reference Standard & Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25×), Substrate Reagent, Stop Solution, Plate Sealer, Product Description, Certificate of Analysis |
| <b>Standard Curve Range</b>    | 0.78 - 50 ng/mL                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sensitivity</b>             | 0.47 ng/mL                                                                                                                                                                                                                                                                                                                                                         |
| <b>Inter Assay</b>             | CV% < 4.5%                                                                                                                                                                                                                                                                                                                                                         |
| <b>Intra Assay</b>             | CV% < 5.69%                                                                                                                                                                                                                                                                                                                                                        |
| <b>Assay Type</b>              | Sandwich-ELISA                                                                                                                                                                                                                                                                                                                                                     |
| <b>Spiking Recovery</b>        | 85-111%                                                                                                                                                                                                                                                                                                                                                            |
| <b>Suitable Sample Type</b>    | Serum, plasma and other biological fluids                                                                                                                                                                                                                                                                                                                          |
| <b>Sample Volume</b>           | 100 uL                                                                                                                                                                                                                                                                                                                                                             |

**Product Application Details**

|                              |       |
|------------------------------|-------|
| <b>Applications</b>          | ELISA |
| <b>Recommended Dilutions</b> | ELISA |

**Images**

Human COL9A3 ELISA Kit (Colorimetric) - Standard Curve Reference





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-75887**

---

|                    |                                               |
|--------------------|-----------------------------------------------|
| H00001299-Q01-10ug | Recombinant Human COL9A3 GST (N-Term) Protein |
| 210-TA-005         | TNF-alpha [Unconjugated]                      |
| H00001299-D01P     | COL9A3 Antibody                               |
| D6050              | IL-6 [HRP]                                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-75887](http://www.novusbio.com/reviews/submit/NBP2-75887)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

