

# Product Datasheet

## **GOLPH3 Antibody** **NBP2-76352**

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-76352](http://www.novusbio.com/NBP2-76352)

Updated 10/10/2022 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-76352](http://www.novusbio.com/reviews/destination/NBP2-76352)




**NBP2-76352****GOLPH3 Antibody**

<b>Product Information</b>	
<b>Unit Size</b>	100 ul
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 2 - 8 C / 1 year from date of receipt
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.09% Sodium Azide
<b>Isotype</b>	IgG
<b>Purity</b>	Antigen Affinity-purified
<b>Buffer</b>	Tris-buffered Saline containing 0.1% BSA
<b>Product Description</b>	
<b>Host</b>	Rabbit
<b>Gene ID</b>	64083
<b>Gene Symbol</b>	GOLPH3
<b>Species</b>	Human, Mouse
<b>Reactivity Notes</b>	Based on 100% sequence identity, this antibody is predicted to react with Rat
<b>Immunogen</b>	Reacts with immunogen between amino acids 1 and 50
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunoprecipitation
<b>Recommended Dilutions</b>	Western Blot 1:1000, Immunoprecipitation 50-100 ul/mg lysate




## Images

**Western Blot: GOLPH3 Antibody [NBP2-76352]** - Western Blot: GOLPH3 Antibody - Samples: Whole cell lysate (15 ug) from HeLa, MCF-7, U2OS, mouse TCMK-1, and mouse NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-GOLPH3 antibody used for WB at 1:1000. Detection: Chemiluminescence with an exposure time of 30 seconds.



**Immunoprecipitation: GOLPH3 Antibody [NBP2-76352]** - Detection of human GOLPH3 by western blot of immunoprecipitates. Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-GOLPH3 antibody NBP2-76352 used for IP at 100 ul per reaction. GOLPH3 was also immunoprecipitated by rabbit anti-GOLPH3 antibody NBP2-76817. For blotting immunoprecipitated GOLPH3, NBP2-76352 was used at 1:1000. Detection: Chemiluminescence with an exposure time of 3 seconds.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-76352](http://www.novusbio.com/reviews/submit/NBP2-76352)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

