

Product Datasheet

YIF1B Antibody NBP2-76381

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-76381

Updated 10/10/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-76381



NBP2-76381**YIF1B Antibody****Product Information**

| | |
|----------------------|---|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 2 - 8 C / 1 year from date of receipt |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Antigen Affinity-purified |
| Buffer | Tris-buffered Saline, 0.1% BSA |

Product Description


| | |
|--------------------|---|
| Host | Rabbit |
| Gene ID | 90522 |
| Gene Symbol | YIF1B |
| Species | Human |
| Immunogen | Immunogen sequence between amino acids 1 and 50 |

Product Application Details


| | |
|------------------------------|--|
| Applications | Western Blot, Immunoprecipitation |
| Recommended Dilutions | Western Blot 1:100, Immunoprecipitation 20-40 ul/mg lysate |

Images

Western Blot: YIF1B Antibody [NBP2-76381] - Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from K-562 cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-YIF1B antibody NBP2-76381 used for IP at 40 ul per reaction. YIF1B was also immunoprecipitated by rabbit anti-YIF1B antibody NBP2-76382. For blotting immunoprecipitated YIF1B, NBP2-76381 was used at 1:100. Detection: Chemiluminescence with an exposure time of 3 seconds.



Western Blot: YIF1B Antibody [NBP2-76381] - Samples: Whole cell lysate (50 ug) from K-562 and A-549 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-YIF1B antibody NBP2-76381 used for WB at 1:1000. Detection: Chemiluminescence with an exposure time of 3 minutes.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-76381

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

