

# Product Datasheet

## **RIPK3/RIP3 Antibody (1312D) [DyLight 680] NBP2-76799FR**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-76799FR](http://www.novusbio.com/NBP2-76799FR)

Updated 9/19/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-76799FR](http://www.novusbio.com/reviews/destination/NBP2-76799FR)



**NBP2-76799FR**

RIPK3/RIP3 Antibody (1312D) [DyLight 680]

| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 0.1 ml                                                                                      |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C in the dark.                                                                    |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | 1312D                                                                                       |
| Preservative            | 0.05% Sodium Azide                                                                          |
| Isotype                 | IgG                                                                                         |
| Conjugate               | DyLight 680                                                                                 |
| Purity                  | Protein A or G purified                                                                     |
| Buffer                  | 50mM Sodium Borate                                                                          |
| Target Molecular Weight | 56.7 kDa                                                                                    |

| Product Description |                                                                                                                                                                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                                                                                                                                                           |
| Gene ID             | 11035                                                                                                                                                                                                                                            |
| Gene Symbol         | RIPK3                                                                                                                                                                                                                                            |
| Species             | Human                                                                                                                                                                                                                                            |
| Immunogen           | RIPK3/RIP3 Antibody (1312D) [DyLight 680] was made to a partial synthetic peptide made to an internal portion of the human RIPK3/RIP3 protein (between amino acids 470-540). This sequence has 100% homology to isoform CRA_a and isoform CRA_b. |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence                   |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/Immunofluorescence                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-76799FR**

---

|                    |                                                   |
|--------------------|---------------------------------------------------|
| NBP2-24891FR       | Rabbit IgG Isotype Control [DyLight 680]          |
| H00011035-P01-25ug | Recombinant Human RIPK3/RIP3 GST (N-Term) Protein |
| 210-TA-005         | TNF-alpha [Unconjugated]                          |
| NBP2-24588PEP      | RIPK3/RIP3 Antibody Blocking Peptide              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-76799FR](http://www.novusbio.com/reviews/submit/NBP2-76799FR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

