

Product Datasheet

VECTASTAIN(R) Universal Quick anti-mouse/rabbit/goat IgG HRP Immunodetection Kit PK-8800-NB

Unit Size: 1 Kit

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/PK-8800-NB

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/PK-8800-NB



PK-8800-NB**VECTASTAIN(R) Universal Quick anti-mouse/rabbit/goat IgG HRP Immunodetection Kit**


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.

Product Description	
Description	<p>Features: Rapid protocol: Staining in less than 20 minutes following primary antibody incubation Working solutions can be used immediately after dilution. High sensitivity, low background. Biotinylated, Universal Pan-Specific secondary antibody recognizes mouse, rabbit, and goat primary antibodies as well as those from related species such as rat, bovine, and sheep.* Available in Ready-To-Use (R.T.U) format: Prediluted, stabilized working solutions of VECTASTAIN® Universal Quick Kit reagents provide the same high sensitivity and low background as the concentrated VECTASTAIN® Universal Kit.</p> <p>*Do not use the pan-specific secondary antibody to stain rat, mouse or other rodent, rabbit, goat, bovine, or sheep tissue due to potential reactivity with endogenous IgG.</p>
Species	Mouse, Goat, Rabbit
Kit Components	1.2 ml peroxidase-streptavidin complex, 2.2 ml biotinylated universal pan-specific antibody (anti-mouse/rabbit/goat IgG), 6 ml normal horse serum
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.

Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry, Immunocytochemistry/ Immunofluorescence
Application Notes	Sufficient reagents to stain approximately 250-500 sections.

Images

Immunohistochemistry: VECTASTAIN(R) Universal Quick anti-mouse/rabbit/goat IgG HRP Immunodetection Kit [PK-8800-NB] - Human Breast Carcinoma: Ki67 (m), VECTASTAIN Universal Quick Kit with DAB (brown) substrate. Hematoxylin (blue) counterstain.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/PK-8800-NB

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

