

Product Datasheet

CaptSure DIY ELISA Kit 6300020

Unit Size: 5 x 96 Tests

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/6300020

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/6300020



6300020**CaptSure DIY ELISA Kit**

Product Information	
Unit Size	5 x 96 Tests
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Description	<p>CaptSure DIY ELISA is the latest cutting-edge ELISA kit with a built-in system for antibody labeling, enabling the rapid and easy development of in-house ELISAs using any antibody pair, while eliminating the need to perform plate coating procedures.</p> <p>The CaptSure DIY ELISA kit has a ready-to-go CaptSure assay plate, pre-coated with our CaptSure antibody, which is specific for the CaptSure peptide. The capture antibody is quickly and easily conjugated with the Lightning-Link CaptSure peptide supplied within the kit. Once added to the plate, the interaction between the CaptSure antibody and the CaptSure peptide immobilizes the immune-complex to the plate.</p> <p>The detection antibody is similarly conjugated with HRP using the Lightning-Link HRP conjugation reagent, also supplied in the kit. Thus, with minimal hands-on time, you can set up both of your antibodies to ensure the best performance of your ELISA.</p>
Kit Components	HRP-Labelled Antibody Dilution Buffer, CaptSure Peptide tagged Antibody Dilution Buffer, Detection Reagent, Stop Solution, Wash Buffer, Sample Dilution Buffer (Secreted Proteins) - 2 EA, Lysis Extraction Buffer (5X) (Cellular Proteins), CaptSure Assay Plates - 5 EA, Lightning-Link CaptSure Peptide Conjugation Kit - 1 EA, Lightning-Link HRP Conjugation Kit - 1 EA
Notes	<p>This product is manufactured by Abcam and distributed by Novus Biologicals.</p> <p>This product is for research use only and is not approved for use in humans or in clinical diagnosis. This product is guaranteed for 1 year from date of receipt and this statement overrides any mentioned guarantee period on the limitations section of this products datasheet. Please contact technical@novusbio.com with questions.</p>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/6300020

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

