

Product Datasheet

p38 gamma/SAPK3 Antibody [Allophycocyanin] AF1644A

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF1644A

Updated 5/14/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF1644A



AF1644A

p38 gamma/SAPK3 Antibody [Allophycocyanin]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	Allophycocyanin
Purity	Antigen Affinity-purified
Buffer	PBS
Product Description	
Host	Goat
Gene ID	6300
Gene Symbol	MAPK12
Species	Human, Mouse, Rat
Specificity/Sensitivity	Detects human, mouse, and rat p38 γ in Western blots. The antibody is highly specific for p38 γ and in Western blots, shows extremely weak cross-reactivity with recombinant p38 α , p38 β , or p38 δ after prolonged exposure to film.
Immunogen	<i>E. coli</i> -derived recombinant human p38 γ Met1-Leu367 Accession # P53778
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF1644A

IC108A	Goat IgG Isotype Control [Allophycocyanin]
NBP1-72551-100ug	Recombinant Human p38 gamma/SAPK3 His Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP2-04236	p38 gamma/SAPK3 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF1644A

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

