

# Product Datasheet

## IL-5 Antibody [DyLight 550] AF2470L

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF2470L](http://www.novusbio.com/AF2470L)

Updated 7/20/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF2470L](http://www.novusbio.com/reviews/destination/AF2470L)



**AF2470L**

IL-5 Antibody [DyLight 550]

| Product Information         |                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.       |
| Storage                     | Store at 4C in the dark.                                                                          |
| Clonality                   | Polyclonal                                                                                        |
| Preservative                | 0.05% Sodium Azide                                                                                |
| Isotype                     | IgG                                                                                               |
| Conjugate                   | DyLight 550                                                                                       |
| Purity                      | Antigen Affinity-purified                                                                         |
| Buffer                      | 50mM Sodium Borate                                                                                |
| Product Description         |                                                                                                   |
| Host                        | Goat                                                                                              |
| Gene ID                     | 3567                                                                                              |
| Gene Symbol                 | IL5                                                                                               |
| Species                     | Equine                                                                                            |
| Specificity/Sensitivity     | Detects equine IL-5 in direct ELISAs and Western blots.                                           |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant equine IL-5<br>Leu20-Gly134<br>Accession # O02699 |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                 |
| Product Application Details |                                                                                                   |
| Applications                | Western Blot, ELISA Detection (Matched Antibody Pair), Immunocytochemistry, Neutralization        |
| Recommended Dilutions       | Western Blot, Neutralization, ELISA Detection (Matched Antibody Pair), Immunocytochemistry        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF2470L**

---

|                |                                        |
|----------------|----------------------------------------|
| AB-108-CL      | Goat IgG Isotype Control [DyLight 550] |
| NBP2-35132-5ug | Recombinant Mouse IL-5 Protein         |
| 210-TA-005     | TNF-alpha [Unconjugated]               |
| M5000          | IL-5 [HRP]                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF2470L](http://www.novusbio.com/reviews/submit/AF2470L)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

