

# Product Datasheet

## MyD88 Antibody [DyLight 350] AF3109D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF3109D](http://www.novusbio.com/AF3109D)

Updated 7/20/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF3109D](http://www.novusbio.com/reviews/destination/AF3109D)



**AF3109D**

MyD88 Antibody [DyLight 350]

| Product Information         |                                                                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                             |
| Storage                     | Store at 4C in the dark.                                                                                                                |
| Clonality                   | Polyclonal                                                                                                                              |
| Preservative                | 0.05% Sodium Azide                                                                                                                      |
| Isotype                     | IgG                                                                                                                                     |
| Conjugate                   | DyLight 350                                                                                                                             |
| Purity                      | Antigen Affinity-purified                                                                                                               |
| Buffer                      | 50mM Sodium Borate                                                                                                                      |
| Product Description         |                                                                                                                                         |
| Host                        | Goat                                                                                                                                    |
| Gene ID                     | 4615                                                                                                                                    |
| Gene Symbol                 | MYD88                                                                                                                                   |
| Species                     | Mouse, Rat                                                                                                                              |
| Specificity/Sensitivity     | Detects mouse and rat MyD88 in Western blots.                                                                                           |
| Immunogen                   | <i>E. coli</i> -derived recombinant mouse MyD88<br>Met1-Pro296<br>Accession # Q3U7M4                                                    |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                       |
| Product Application Details |                                                                                                                                         |
| Applications                | Western Blot, Immunoprecipitation, Proximity Ligation Assay, CyTOF-ready, Immunocytochemistry, Intracellular Staining by Flow Cytometry |
| Recommended Dilutions       | Western Blot, Immunoprecipitation, Proximity Ligation Assay, Intracellular Staining by Flow Cytometry, Immunocytochemistry, CyTOF-ready |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF3109D**

---

|                |                                        |
|----------------|----------------------------------------|
| AB-108-CD      | Goat IgG Isotype Control [DyLight 350] |
| NB100-56698PEP | MyD88 Antibody Blocking Peptide        |
| 210-TA-005     | TNF-alpha [Unconjugated]               |
| NBP2-29328     | MyD88 Inhibitor Peptide Set            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF3109D](http://www.novusbio.com/reviews/submit/AF3109D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

