

# Product Datasheet

## M-Ras/R-Ras3 Antibody [FITC] AF3606F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF3606F](http://www.novusbio.com/AF3606F)

Updated 5/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF3606F](http://www.novusbio.com/reviews/destination/AF3606F)



**AF3606F****M-Ras/R-Ras3 Antibody [FITC]**

| Product Information         |                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                    |
| Storage                     | Store at 4C in the dark.                                                                                       |
| Clonality                   | Polyclonal                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                             |
| Isotype                     | IgG                                                                                                            |
| Conjugate                   | FITC                                                                                                           |
| Purity                      | Antigen Affinity-purified                                                                                      |
| Buffer                      | PBS                                                                                                            |
| Product Description         |                                                                                                                |
| Host                        | Goat                                                                                                           |
| Gene ID                     | 22808                                                                                                          |
| Gene Symbol                 | MRAS                                                                                                           |
| Species                     | Human, Rat                                                                                                     |
| Specificity/Sensitivity     | Detects human and rat M-Ras/R-Ras3 in Western blots. Does not detect recombinant human R-Ras or TC21 (R-Ras2). |
| Immunogen                   | <i>E. coli</i> -derived recombinant human M-Ras Met1-Leu208 (Lys181Glu)<br>Accession # O14807                  |
| Product Application Details |                                                                                                                |
| Applications                | Western Blot                                                                                                   |
| Recommended Dilutions       | Western Blot                                                                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF3606F**

---

|                    |                                            |
|--------------------|--------------------------------------------|
| NBP2-33376H        | Blue Marker Antibody (6F4-F6) [HRP]        |
| NB410-28088F-500ul | Goat IgG Isotype Control [FITC]            |
| NBP1-78834-100ug   | Recombinant Human M-Ras/R-Ras3 His Protein |
| 236-EG-200         | EGF [Unconjugated]                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF3606F](http://www.novusbio.com/reviews/submit/AF3606F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

