

# Product Datasheet

## Fas/TNFRSF6/CD95 Antibody [DyLight 755] AF435Z

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF435Z](http://www.novusbio.com/AF435Z)

Updated 7/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF435Z](http://www.novusbio.com/reviews/destination/AF435Z)



**AF435Z**

Fas/TNFRSF6/CD95 Antibody [DyLight 755]

| Product Information         |                                                                                                                                                                                |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                       |
| Clonality                   | Polyclonal                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                             |
| Isotype                     | IgG                                                                                                                                                                            |
| Conjugate                   | DyLight 755                                                                                                                                                                    |
| Purity                      | Antigen Affinity-purified                                                                                                                                                      |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                |
| Host                        | Goat                                                                                                                                                                           |
| Gene ID                     | 355                                                                                                                                                                            |
| Gene Symbol                 | FAS                                                                                                                                                                            |
| Species                     | Mouse                                                                                                                                                                          |
| Specificity/Sensitivity     | Detects mouse Fas/TNFRSF6/CD95 in ELISAs and Western blots.                                                                                                                    |
| Immunogen                   | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant mouse Fas/TNFRSF6/CD95<br>Gly14-Arg169<br>Accession # P25446                                           |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                              |
| Product Application Details |                                                                                                                                                                                |
| Applications                | Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair) |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry, ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), CyTOF-ready |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF435Z**

---

|                |                                            |
|----------------|--------------------------------------------|
| AB-108-CZ      | Goat IgG Isotype Control [DyLight 755]     |
| NBP2-61594-5ug | Recombinant Human Fas/TNFRSF6/CD95 Protein |
| 210-TA-005     | TNF-alpha [Unconjugated]                   |
| 7398-FS-050    | Fas/TNFRSF6/CD95                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF435Z](http://www.novusbio.com/reviews/submit/AF435Z)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

