

# Product Datasheet

## PLD1 Antibody [HRP] AF5615H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF5615H](http://www.novusbio.com/AF5615H)

Updated 5/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF5615H](http://www.novusbio.com/reviews/destination/AF5615H)



**AF5615H**

PLD1 Antibody [HRP]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Polyclonal                                                                                  |
| Preservative                | No Preservative                                                                             |
| Isotype                     | IgG                                                                                         |
| Conjugate                   | HRP                                                                                         |
| Purity                      | Antigen Affinity-purified                                                                   |
| Buffer                      | PBS                                                                                         |
| Product Description         |                                                                                             |
| Host                        | Goat                                                                                        |
| Gene ID                     | 5337                                                                                        |
| Gene Symbol                 | PLD1                                                                                        |
| Species                     | Human                                                                                       |
| Specificity/Sensitivity     | Detects endogenous human PLD1 in Western blots.                                             |
| Immunogen                   | <i>E. coli</i> -derived recombinant human PLD1<br>Met1-Pro140<br>Accession # Q13393         |
| Product Application Details |                                                                                             |
| Applications                | Western Blot                                                                                |
| Recommended Dilutions       | Western Blot                                                                                |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF5615H**

---

|                    |                                             |
|--------------------|---------------------------------------------|
| AB-108-CH          | Goat IgG Isotype Control [HRP]              |
| H00005337-Q01-10ug | Recombinant Human PLD1 GST (N-Term) Protein |
| H00005555-P01-10ug | Recombinant Human PRH2 GST (N-Term) Protein |
| NBP3-39372         | Human PLD1 ELISA Kit (Colorimetric)         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF5615H](http://www.novusbio.com/reviews/submit/AF5615H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

