

# Product Datasheet

## Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 Antibody [DyLight 350] AF6664D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF6664D](http://www.novusbio.com/AF6664D)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF6664D](http://www.novusbio.com/reviews/destination/AF6664D)



**AF6664D**

Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 Antibody [DyLight 350]

| Product Information         |                                                                                                                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                         |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                            |
| Clonality                   | Polyclonal                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                  |
| Isotype                     | IgG                                                                                                                                                                                                 |
| Conjugate                   | DyLight 350                                                                                                                                                                                         |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                           |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                  |
| Product Description         |                                                                                                                                                                                                     |
| Host                        | Sheep                                                                                                                                                                                               |
| Gene ID                     | 11041                                                                                                                                                                                               |
| Gene Symbol                 | B4GAT1                                                                                                                                                                                              |
| Species                     | Human, Mouse                                                                                                                                                                                        |
| Specificity/Sensitivity     | Detects human and mouse $\beta$ -1,4-Glucuronyltransferase 1/B4GAT1 in Western blot. In direct ELISAs, less than 1% cross-reactivity with recombinant human (rh) B4GalT1 and rhB3GNT 2 is observed. |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human $\beta$ -1,4-Glucuronyltransferase 1/B4GAT1<br>Asp43-Cys415<br>Accession # O43505                                                     |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                   |
| Product Application Details |                                                                                                                                                                                                     |
| Applications                | Western Blot                                                                                                                                                                                        |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                              |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF6664D**

---

|               |                                                                               |
|---------------|-------------------------------------------------------------------------------|
| NBP1-88952PEP | Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 Recombinant Protein Antigen |
| MEP00B        | Erythropoietin/EPO [HRP]                                                      |
| 6664-GT-050   | Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 [Unconjugated]              |
| NB100-524     | NOD2 Antibody (2D9) - BSA Free                                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF6664D](http://www.novusbio.com/reviews/submit/AF6664D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

