

# Product Datasheet

## TLE2 Antibody [Janelia Fluor® 549] AF7216I

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF7216I](http://www.novusbio.com/AF7216I)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF7216I](http://www.novusbio.com/reviews/destination/AF7216I)



**AF7216I**

TLE2 Antibody [Janelia Fluor® 549]

| Product Information         |                                                                                                                                                        |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                               |
| Clonality                   | Polyclonal                                                                                                                                             |
| Preservative                | 0.05% Sodium Azide                                                                                                                                     |
| Isotype                     | IgG                                                                                                                                                    |
| Conjugate                   | Janelia Fluor 549                                                                                                                                      |
| Purity                      | Antigen Affinity-purified                                                                                                                              |
| Buffer                      | 50mM Sodium Borate                                                                                                                                     |
| Product Description         |                                                                                                                                                        |
| Host                        | Sheep                                                                                                                                                  |
| Gene ID                     | 7089                                                                                                                                                   |
| Gene Symbol                 | TLE2                                                                                                                                                   |
| Species                     | Human                                                                                                                                                  |
| Specificity/Sensitivity     | Detects human TLE2 in direct ELISAs. In direct ELISAs, less than 2% cross-reactivity with recombinant human (rh) TLE1, rhTLE3, and rhTLE4 is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant human TLE2<br>Ile11-Ser193<br>Accession # Q04725                                                                   |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                  |
| Product Application Details |                                                                                                                                                        |
| Applications                | Immunohistochemistry                                                                                                                                   |
| Recommended Dilutions       | Immunohistochemistry                                                                                                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                 |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to AF7216I**

---

|                    |                                             |
|--------------------|---------------------------------------------|
| H00007089-P01-10ug | Recombinant Human TLE2 GST (N-Term) Protein |
| M6000B-1           | IL-6 [HRP]                                  |
| NBL1-16945         | TLE2 Overexpression Lysate                  |
| AF808              | Osteopontin/OPN Antibody [Unconjugated]     |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF7216I](http://www.novusbio.com/reviews/submit/AF7216I)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

