

Product Datasheet

PI 3-Kinase C2 beta Antibody [DyLight 550] AF7249L

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF7249L

Updated 5/28/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF7249L



AF7249L

PI 3-Kinase C2 beta Antibody [DyLight 550]

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 550
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate

Product Description

Host	Goat
Gene ID	5287
Gene Symbol	PIK3C2B
Species	Human
Specificity/Sensitivity	Detects human PI 3□Kinase C2β in direct ELISAs and Western blots.
Immunogen	<i>E. coli</i> -derived recombinant human PI 3□Kinase C2β Ser2-Ser134 Accession # O00750
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details

Applications	Western Blot, Immunocytochemistry
Recommended Dilutions	Western Blot, Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF7249L

NB410-28088R	Goat IgG Isotype Control [DyLight 550]
NBP2-58352PEP	PI 3-Kinase C2 beta Recombinant Protein Antigen
236-EG-200	EGF [Unconjugated]
NBP1-86795PEP	PI 3-Kinase C2 beta Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF7249L

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

