

# Product Datasheet

## 4-1BB/TNFRSF9/CD137 Antibody [Janelia Fluor® 549] AF838I

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF838I](http://www.novusbio.com/AF838I)

Updated 7/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF838I](http://www.novusbio.com/reviews/destination/AF838I)



**AF838I****4-1BB/TNFRSF9/CD137 Antibody [Janelia Fluor® 549]**

| <b>Product Information</b>         |                                                                                                                                                                                                                       |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                           |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                              |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                                            |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                    |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                   |
| <b>Conjugate</b>                   | Janelia Fluor 549                                                                                                                                                                                                     |
| <b>Purity</b>                      | Antigen Affinity-purified                                                                                                                                                                                             |
| <b>Buffer</b>                      | 50mM Sodium Borate                                                                                                                                                                                                    |
| <b>Product Description</b>         |                                                                                                                                                                                                                       |
| <b>Host</b>                        | Goat                                                                                                                                                                                                                  |
| <b>Gene ID</b>                     | 3604                                                                                                                                                                                                                  |
| <b>Gene Symbol</b>                 | TNFRSF9                                                                                                                                                                                                               |
| <b>Species</b>                     | Human                                                                                                                                                                                                                 |
| <b>Specificity/Sensitivity</b>     | Detects 4-1BB/TNFRSF9/CD137 in direct ELISAs and Western blots.                                                                                                                                                       |
| <b>Immunogen</b>                   | Mouse myeloma cell line NS0-derived recombinant human 4-1BB/TNFRSF9/CD137 Leu24-Gln186 Accession # Q07011                                                                                                             |
| <b>Notes</b>                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                                                                 |
| <b>Product Application Details</b> |                                                                                                                                                                                                                       |
| <b>Applications</b>                | Western Blot, Flow Cytometry, Immunohistochemistry, Agonist Activity, CyTOF-ready, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair), Immunocytochemistry |
| <b>Recommended Dilutions</b>       | Western Blot, Flow Cytometry, Immunohistochemistry, Agonist Activity, ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), Immunocytochemistry, CyTOF-ready |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF838I**

---

|                |                                               |
|----------------|-----------------------------------------------|
| AB-108-CI      | Goat IgG Isotype Control [Janelia Fluor 549]  |
| NBP2-34958-5ug | Recombinant Human 4-1BB/TNFRSF9/CD137 Protein |
| 210-TA-005     | TNF-alpha [Unconjugated]                      |
| 937-4B-050     | 4-1BB/TNFRSF9/CD137                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF838I](http://www.novusbio.com/reviews/submit/AF838I)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

