

# Product Datasheet

## SIGIRR Antibody (161917) [PE] FAB1092P

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1092P](http://www.novusbio.com/FAB1092P)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1092P](http://www.novusbio.com/reviews/destination/FAB1092P)



**FAB1092P**

SIGIRR Antibody (161917) [PE]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.5 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 161917                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | PE                                                                                          |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                              |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rat                                                                                                                                                                                                                                                                                          |
| Gene ID                 | 59307                                                                                                                                                                                                                                                                                        |
| Gene Symbol             | SIGIRR                                                                                                                                                                                                                                                                                       |
| Species                 | Mouse                                                                                                                                                                                                                                                                                        |
| Specificity/Sensitivity | Detects mouse SIGIRR in direct ELISAs and Western blots. In direct ELISAs, this antibody does not cross-react with recombinant mouse (rm) IL-1 R1, rmIL-1 R2, rhIL-1 RAcP, recombinant human (rh) IL-1 R4, rmIL-18 R, recombinant rat IL-1 R6, rmIL-1 R7, rhIL-1 R8, rhIL-1 R9, or rhSIGIRR. |
| Immunogen               | Mouse myeloma cell line NS0-derived recombinant mouse SIGIRR Met1-His117<br>Accession # Q9JLZ8                                                                                                                                                                                               |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB1092P**

---

|                |                                  |
|----------------|----------------------------------|
| NB100-56890PEP | SIGIRR Antibody Blocking Peptide |
| 210-TA-005     | TNF-alpha [Unconjugated]         |
| 992-SG-100     | SIGIRR [Unconjugated]            |
| M6000B-1       | IL-6 [HRP]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1092P](http://www.novusbio.com/reviews/submit/FAB1092P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

