

Product Datasheet

Siglec-3/CD33 Antibody (2337A) [Allophycocyanin] FAB11373A

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB11373A

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB11373A



FAB11373A

Siglec-3/CD33 Antibody (2337A) [Allophycocyanin]

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	2337A
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	Allophycocyanin
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS

Product Description

Host	Rabbit
Gene ID	945
Species	Human
Specificity/Sensitivity	Detects human Siglec-3/CD33 in direct ELISAs.
Immunogen	Mouse myeloma cell line, NS0-derived human Siglec-3/CD33 Asp18-His259 Accession # P20138

Product Application Details

Applications	Flow Cytometry, CyTOF-ready
Recommended Dilutions	Flow Cytometry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

Siglec-3/CD33 Antibody (2337A) [Allophycocyanin] [FAB11373A] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB11373A

IC1051A	Rabbit IgG Isotype Control (60024B) [Allophycocyanin]
NBP2-32709PEP	Siglec-3/CD33 Recombinant Protein Antigen
DC140	CD14 [HRP]
1137-SL-050	Siglec-3/CD33 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB11373A

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

