

# Product Datasheet

## **CGRP2 Antibody (1057805) [DyLight 350] FAB11382D**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB11382D](http://www.novusbio.com/FAB11382D)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB11382D](http://www.novusbio.com/reviews/destination/FAB11382D)



**FAB11382D**

CGRP2 Antibody (1057805) [DyLight 350]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 1057805                                                                                     |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2b                                                                                       |
| Conjugate           | DyLight 350                                                                                 |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                              |
|-------------------------|----------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                        |
| Gene ID                 | 797                                                                                          |
| Species                 | Human                                                                                        |
| Specificity/Sensitivity | Detects Human Procalcitonin in direct ELISA.                                                 |
| Immunogen               | Human CGRP synthetic peptide (Tocris (Catalog # <a href="#">3012</a> )<br>Accession # P10092 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.            |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, ELISA, Immunocytochemistry                               |
| Recommended Dilutions       | Western Blot, ELISA, Immunocytochemistry                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB11382D**

---

|              |                                                    |
|--------------|----------------------------------------------------|
| NBP2-27231UV | Mouse IgG2b Isotype Control (MPC-11) [DyLight 350] |
| NBP1-30052   | Choline Acetyltransferase/ChAT Antibody            |
| MAB1417      | Insulin Antibody (182410) [Unconjugated]           |
| NBP1-46535   | Neuropeptide Y Antibody                            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB11382D](http://www.novusbio.com/reviews/submit/FAB11382D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

