

# Product Datasheet

## Serpin C1/Antithrombin-III Antibody (208730) [HRP] FAB1287H

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1287H](http://www.novusbio.com/FAB1287H)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1287H](http://www.novusbio.com/reviews/destination/FAB1287H)



**FAB1287H**

Serpine C1/Antithrombin-III Antibody (208730) [HRP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.5 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 208730                                                                                      |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG2a                                                                                       |
| Conjugate           | HRP                                                                                         |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                            |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rat                                                                                                                                                                                                                                                                                        |
| Gene ID                 | 462                                                                                                                                                                                                                                                                                        |
| Gene Symbol             | SERPINC1                                                                                                                                                                                                                                                                                   |
| Species                 | Mouse                                                                                                                                                                                                                                                                                      |
| Specificity/Sensitivity | Detects mouse Serpin C1/Antithrombin-III in direct ELISAs and Western blots. In direct ELISAs and Western blots, 100% cross-reactivity with recombinant human (rh) Serpin C1 is observed and no cross-reactivity with rhSerpine A1, A3, F2, recombinant mouse Serpin D1 or A5 is observed. |
| Immunogen               | Mouse myeloma cell line NS0-derived recombinant mouse Serpin C1/Antithrombin-III<br>His33-Asn465<br>Accession # P32261                                                                                                                                                                     |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB1287H**

---

|               |                                                        |
|---------------|--------------------------------------------------------|
| NBP1-90076PEP | Serpin C1/Antithrombin-III Recombinant Protein Antigen |
| DCF300        | Coagulation Factor III/Tissue Factor [HRP]             |
| 1267-PI-010   | Serpin C1/Antithrombin-III [Unconjugated]              |
| MAB1455       | Albumin Antibody (188835) [Unconjugated] - Serum       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1287H](http://www.novusbio.com/reviews/submit/FAB1287H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

