

# Product Datasheet

## Kininogen Antibody (236006) [DyLight 350] FAB15691D

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB15691D](http://www.novusbio.com/FAB15691D)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB15691D](http://www.novusbio.com/reviews/destination/FAB15691D)



**FAB15691D**

Kininogen Antibody (236006) [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Clone                       | 236006                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Gene ID                     | 3827                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Gene Symbol                 | KNG1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Specificity/Sensitivity     | Detects human Kininogen Light Chain in direct ELISAs and Western blots. In Western blots, it reacts with the full length (amino acids 19–644) and two forms of the light chain starting at amino acids 390 and 438, respectively, and does not react with the <i>E. coli</i> -expressed D5/Kininostatin. In Western blots, no cross-reactivity with recombinant human (rh) Cystatins A, B, C, D, E/M, F, S, SA, SN, rhFetuin A, B, rhHPRG, or recombinant mouse Kininogen is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Kininogen Gln19-Ser644 (Ile581Thr) (predicted)<br>Accession # P01042.2                                                                                                                                                                                                                                                                                                                                                           |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                                                                                      |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB15691D**

---

|               |                                                    |
|---------------|----------------------------------------------------|
| NBP2-27231UV  | Mouse IgG2b Isotype Control (MPC-11) [DyLight 350] |
| NBP1-89993PEP | Kininogen Recombinant Protein Antigen              |
| 210-TA-005    | TNF-alpha [Unconjugated]                           |
| 1569-PI-010   | Kininogen                                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB15691D](http://www.novusbio.com/reviews/submit/FAB15691D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

