

# Product Datasheet

## Carbonic Anhydrase X/CA10 Antibody (862319) [Biotin] FAB2189B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB2189B](http://www.novusbio.com/FAB2189B)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB2189B](http://www.novusbio.com/reviews/destination/FAB2189B)



**FAB2189B**

Carbonic Anhydrase X/CA10 Antibody (862319) [Biotin]

| Product Information         |                                                                                                                      |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                               |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                          |
| Storage                     | Store at 4C in the dark.                                                                                             |
| Clonality                   | Monoclonal                                                                                                           |
| Clone                       | 862319                                                                                                               |
| Preservative                | 0.05% Sodium Azide                                                                                                   |
| Isotype                     | IgG2b                                                                                                                |
| Conjugate                   | Biotin                                                                                                               |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                           |
| Buffer                      | PBS                                                                                                                  |
| Product Description         |                                                                                                                      |
| Host                        | Mouse                                                                                                                |
| Gene ID                     | 56934                                                                                                                |
| Gene Symbol                 | CA10                                                                                                                 |
| Species                     | Human                                                                                                                |
| Specificity/Sensitivity     | Detects human Carbonic Anhydrase X/CA10 in ELISAs and Western blots.                                                 |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Carbonic Anhydrase X/CA10<br>Met1-Asn300<br>Accession # Q9NS85 |
| Product Application Details |                                                                                                                      |
| Applications                | Western Blot                                                                                                         |
| Recommended Dilutions       | Western Blot                                                                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB2189B**

---

|               |                                                       |
|---------------|-------------------------------------------------------|
| NBP2-27231B   | Mouse IgG2b Isotype Control (MPC-11) [Biotin]         |
| NBP2-31586PEP | Carbonic Anhydrase X/CA10 Recombinant Protein Antigen |
| 2189-CA-010   | Carbonic Anhydrase X/CA10 [Unconjugated]              |
| H00006690-M01 | SPINK1 Antibody (4D4) - Azide and BSA Free            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB2189B](http://www.novusbio.com/reviews/submit/FAB2189B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

