

# Product Datasheet

## Kininogen Antibody (288422) [PerCP] FAB2206C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB2206C](http://www.novusbio.com/FAB2206C)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB2206C](http://www.novusbio.com/reviews/destination/FAB2206C)



**FAB2206C**

Kininogen Antibody (288422) [PerCP]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                              |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                               |
| Clone                       | 288422                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                       |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Conjugate                   | PerCP                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                                                                               |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Host                        | Rat                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Gene ID                     | 3827                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Gene Symbol                 | Kng1                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Species                     | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Specificity/Sensitivity     | Detects mouse Kininogen in direct ELISAs and Western blots. In Western blots, reacts with the full length mouse Kininogen (aa 21-661), but not with the light chain alone. In Western blots, no cross-reactivity with recombinant human (rh) Cystatin D, F, SA, SN, recombinant mouse (rm) Cystatin A, B, C, E/M, rmCystatin/Stefin Homolog, rmFetuin A, rhFetuin B, rmHPRG, rhKininogen, or rhKininostatin is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant mouse Kininogen Glu21-Ser661<br>Accession # O08677                                                                                                                                                                                                                                                                                                                       |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                                                                        |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                                                                        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB2206C**

---

|               |                                                  |
|---------------|--------------------------------------------------|
| NBP1-89994PEP | Kininogen Recombinant Protein Antigen            |
| 210-TA-005    | TNF-alpha [Unconjugated]                         |
| 1569-PI-010   | Kininogen                                        |
| MAB1455       | Albumin Antibody (188835) [Unconjugated] - Serum |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB2206C](http://www.novusbio.com/reviews/submit/FAB2206C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

