

# Product Datasheet

## Serpin D1/Heparin Cofactor II Antibody (373008) [HRP] FAB3198H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB3198H](http://www.novusbio.com/FAB3198H)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB3198H](http://www.novusbio.com/reviews/destination/FAB3198H)



**FAB3198H**

Serpine D1/Heparin Cofactor II Antibody (373008) [HRP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 373008                                                                                      |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG2b                                                                                       |
| Conjugate           | HRP                                                                                         |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                     |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                                                                                                               |
| Gene ID                 | 3053                                                                                                                                                                                                                                                |
| Gene Symbol             | SERPIND1                                                                                                                                                                                                                                            |
| Species                 | Human                                                                                                                                                                                                                                               |
| Specificity/Sensitivity | Detects human Serpin D1/Heparin Cofactor II in direct ELISAs and Western blots. In Western blots, no cross-reactivity with recombinant human Serpin A3, A4, B6, B8, C1, I1, I2, recombinant mouse Serpin A1, A6, D1, E1, E2, F1, or F2 is observed. |
| Immunogen               | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human Serpin D1/Heparin Cofactor II<br>Gly20-Arg498<br>Accession # P05546                                                                                                   |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Immunoprecipitation                                      |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB3198H**

|                    |                                                           |
|--------------------|-----------------------------------------------------------|
| NBP2-27231H        | Mouse IgG2b Isotype Control (MPC-11) [HRP]                |
| NBP2-33764PEP      | Serpin D1/Heparin Cofactor II Recombinant Protein Antigen |
| H00005555-P01-10ug | Recombinant Human PRH2 GST (N-Term) Protein               |
| 3198-PI-010        | Serpin D1/Heparin Cofactor II                             |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB3198H](http://www.novusbio.com/reviews/submit/FAB3198H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

