

Product Datasheet

Methionine Aminopeptidase 1/METAP1 Antibody (405419) **[Janelia Fluor® 669]** **FAB3537JF669**

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB3537JF669

Updated 5/10/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB3537JF669



FAB3537JF669

Methionine Aminopeptidase 1/METAP1 Antibody (405419) [Janelia Fluor® 669]

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 405419 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2b |
| Conjugate | Janelia Fluor 669 |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Mouse |
| Gene ID | 23173 |
| Gene Symbol | METAP1 |
| Species | Human |
| Specificity/Sensitivity | Detects human Methionine Aminopeptidase 1/METAP1 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human Methionine Aminopeptidase 2 is observed. |
| Immunogen | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human Methionine Aminopeptidase 1/METAP1 His52-Phe386 Accession # P53582 |
| Notes | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus. |
| Product Application Details | |
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB3537JF669

| | |
|------------------|--|
| NBP1-78803-100ug | Recombinant E. coli Methionine Aminopeptidase 1/METAP1 His Protein |
| 1290-IL-010 | IL-23A/IL-23 P19 [Unconjugated] |
| 3537-ZN-010 | Methionine Aminopeptidase 1/METAP1 [Unconjugated] |
| NB100-2331 | LC3A Antibody - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB3537JF669

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

