

# Product Datasheet

## **Sulfatase Modifying Factor 2/SUMF2 Antibody (330312)** **[DyLight 405]** **FAB3594E**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB3594E](http://www.novusbio.com/FAB3594E)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB3594E](http://www.novusbio.com/reviews/destination/FAB3594E)



**FAB3594E**

Sulfatase Modifying Factor 2/SUMF2 Antibody (330312) [DyLight 405]

| Product Information         |                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                        |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                   |
| Storage                     | Store at 4C in the dark.                                                                                                                                                      |
| Clonality                   | Monoclonal                                                                                                                                                                    |
| Clone                       | 330312                                                                                                                                                                        |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                            |
| Isotype                     | IgG2b                                                                                                                                                                         |
| Conjugate                   | DyLight 405                                                                                                                                                                   |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                    |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                            |
| Product Description         |                                                                                                                                                                               |
| Host                        | Mouse                                                                                                                                                                         |
| Gene ID                     | 25870                                                                                                                                                                         |
| Gene Symbol                 | SUMF2                                                                                                                                                                         |
| Species                     | Human                                                                                                                                                                         |
| Specificity/Sensitivity     | Detects human Sulfatase Modifying Factor 2/SUMF2 in direct ELISAs and Western blots. No cross-reactivity with recombinant human SUMF1 or recombinant mouse SUMF2 is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Sulfatase Modifying Factor 2/SUMF2<br>Gln26-Leu301<br>Accession # Q8NBJ7                                                |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                             |
| Product Application Details |                                                                                                                                                                               |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                             |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                        |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB3594E**

|               |                                                                |
|---------------|----------------------------------------------------------------|
| NBP2-27231V   | Mouse IgG2b Isotype Control (MPC-11) [DyLight 405]             |
| NBP1-84145PEP | Sulfatase Modifying Factor 2/SUMF2 Recombinant Protein Antigen |
| DY413-05      | IL-13 [Biotin]                                                 |
| AF3264        | Cadherin-13 Antibody [Unconjugated]                            |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB3594E](http://www.novusbio.com/reviews/submit/FAB3594E)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

