

# Product Datasheet

## Integrin beta 4/CD104 Antibody (422325) [HRP] FAB4060H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB4060H](http://www.novusbio.com/FAB4060H)

Updated 8/1/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB4060H](http://www.novusbio.com/reviews/destination/FAB4060H)



**FAB4060H**

Integrin beta 4/CD104 Antibody (422325) [HRP]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	422325
Preservative	No Preservative
Isotype	IgG2b
Conjugate	HRP
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	3691
Gene Symbol	ITGB4
Species	Human
Specificity/Sensitivity	Detects human Integrin $\beta$ 4/CD104 in Western blots. In Western blots, no cross-reactivity with recombinant human Integrin $\beta$ 1, $\beta$ 2, $\beta$ 3, $\beta$ 5, $\beta$ 6, $\beta$ 8, or recombinant mouse Integrin $\beta$ 4 is observed.
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant human Integrin $\beta$ 4/CD104 Asn28-Ser710 Accession # P16144
Product Application Details	
Applications	Western Blot, Flow Cytometry, CyTOF-ready, Immunocytochemistry
Recommended Dilutions	Western Blot, Flow Cytometry, Immunocytochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB4060H**

NBP2-27231H	Mouse IgG2b Isotype Control (MPC-11) [HRP]
NBP2-61461-10ug	Recombinant Human Integrin beta 4/CD104 His Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP2-38297PEP	Integrin beta 4/CD104 Recombinant Protein Antigen

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB4060H](http://www.novusbio.com/reviews/submit/FAB4060H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

